SALVAFAST

ENGLISH



Salvafast upper arm bands *are a fast and easy limb patient restraint intended for post thoracic surgeries. Intended for patients with abdominal wounds, ascites and/or multiple lines which complicate the use of a waist belt.*

Made in comfortable neoprene and double Hook/Loop fastening system that allows a quick application and release. Fixing to the bed by secure closure system to keep the system permanently attached on the bed frame.

ARM BANDS

FEATURES



Safety lock & Magnetic key



Quick attachment Double hook/loop fastening system



Permanent bed installation

SPECIFICATIONS

PRODUCT FAMILY: SALVAFAST

PRODUCT NAME: Arm bands

INTENDED USE: Fast restraint of disturbed, violent, aggressive or disoriented patients

Basic UDI-DI: 843603535SVTP8

SIZE CHART

Ref.	Measures (cm)
SVT3755	45 x 12

COLOUR

Black / Blue

COMPOSITION

80% Polyester - 20% Neoprene
Stainless steel eyelets
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Salvafast arm bands
- 2 - Safety locks
- 1 - Magnetic key
- 1 - Instructions manual

ATTENTION

- Do not use in case of fractures or wounds at the upper arm.
- Always use the upper arm fixation together with wrist restraints.

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE