mobi-tools Glide

ENGLISH



Mobi-tools transfer board *designed to transfer patients from bed to a wheelchair, chair, armchair or vice versa. Significantly reduces the effort required by the caregiver and increases the patient safety due to permanently has a rigid base. Two handles for transportation and handling.*

Made of sturdy and lightweight material, with slick upper face and anti-slip surface on the back.

TRANSFER BOARD 60 x 33 cm


FEATURES

X-RAY / MRI

X-ray translucent and MRI safe



Semi-rigid durable and slick material



Anti-slip tape on back for better grip

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Glide

PRODUCT NAME: Transfer board 60 x 33 cm

INTENDED USE: Transfers of seated patients between two surfaces (bed, chair, couch, wheelchair, toilet, car...)

BASIC UDI-DI: 843603535TMGMQ

SIZE CHART

Ref.	Measures (cm)	Height (mm)
TMG6400	60 x 33	5

COLOUR

- Blue
- Grey tape

COMPOSITION

Polyethylene
LATEX FREE

WARNING

- Maximum weight: 135 Kg.

LAUNDRY INSTRUCTIONS

- Can be cleaned with a cloth dampened with an aqueous solution and neutral detergent. Compatible with surface disinfectant detergents.

PACKAGING

BULK PACKAGING

- 1 - Mobi-tools Glide - Transfer board 60x33

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE

rev. 211229