mobi-tools Glide

ENGLISH



Mobi-tools® Glide transfer board of 75 x 33 cm *is a semi-rigid board for transferring seated patients between two surfaces (bed, chair, armchair, wheelchair, toilet, car, etc.) in a comfortable, easy, fast, safe way and with minimum effort.*

Made from a 6mm thick blue polyethylene board, with a smooth upper surface with sliding properties and two non-slip bands on the lower surface that secure the board on supporting surfaces.

It incorporates two longitudinal slots that, when folded, allow better adaptation to all types of scenarios. Additionally, it has a handle on one side that allows for comfortable transport or hanging for storage.

TRANSFER BOARD 75 x 33 cm

FEATURES

X-RAY / MRI — **X-ray translucent and MRI safe**

 **Semi-rigid durable and slick material**

 **Anti-slip tape on back for better grip**

 **2 bending guides to fit any surface**

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Glide

PRODUCT NAME: Transfer board 75 x 33 cm

INTENDED USE: Transfers of seated patients between two surfaces (bed, chair, couch, wheelchair, toilet, car...)

BASIC UDI-DI: 843603535TMGMQ

SIZE CHART

Ref.	Measures (cm)	Height (mm)	Weight (kg)
TMG6410	75 x 33	6	1,40

COLOUR

- Blue
- Grey tape

COMPOSITION

Polyethylene
LATEX FREE

WARNING

- Maximum weight: 135 KG.

LAUNDRY INSTRUCTIONS

- Can be cleaned with a cloth dampened with an aqueous solution and neutral detergent. Compatible with surface disinfectant detergents.

PACKAGING

BULK PACKAGING

- 1 - Mobi-tools Glide - Transfer board 75x33

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council
*This data sheet could be changed without prior notice

CE

rev. 210518