MEDIROLLER Pro

ENGLISH



Medi-roller PRO transfer board *allows the patient to be moved laterally in a comfortable and safety way, without lifting the patient and with less effort. It consists of a compact, ultralight and highly resistant padded core with a super-slippery surface finish that reduces friction with the tubular sleeve.*

The tubular sleeve has an antistatic surface treatment. It´s easy to sanitize, with a pleasant feel on the outside and slippery surface on the inside, which allows a sliding of the board without the need for great effort.

The material itself does not favor the proliferation of bacteria, non-toxic and hypoallergenic.

Available in 2 versions, large folding board has 2 handles for transport and short one-piece board has 2 handles. The disposable cover for infection prevention is also available for both formats.




TMR6225 - Short board

MEDI-ROLLER PRO

FEATURES



X-ray translucent and MRI safe



European Patent (EP3245996A1)



Reinforced ripstop nylon sliding cover



Thermoformed PE foam core with ultra sliding fabric surface on both sides



Lightweight & handy
Large: 1,4 Kg
Short: 0,70 Kg



2 models:
Large - Foldable - 2 handles
Short - One piece - 2 handles

SPECIFICATIONS

PRODUCT FAMILY: Medi-roller

PRODUCT NAME: Medi-roller PRO

PATENT: EP3245996A1

INTENDED USE: Perform transfers of patients in bed, stretcher, operating table or x-ray...

BASIC UDI-DI: 843603535TMRNE

COEFICIENT OF FRICTION (Static COF): 0,11829586

SIZE CHART

Ref.	Model	Measures (cm)	Weight (Kg)
TMR6205	Foldable	180 x 49	1,40
TMR6225	Short	94 x 49	0,70

COLOUR

- Grey board
- Orange cover

COMPOSITION

Board: Polyethylene
Cover: Polyamide
LATEX FREE

LAUNDRY INSTRUCTIONS

Board: Clean with damp cloth and mild soap or surface disinfectant detergent.

Cover: - Machine wash max. 60ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean



DISINFECTION RECOMMENDATIONS

- For infection prevention and control use the disposable protective cover (ref. TMR6275).
- Use a surface disinfectant or disinfectant wipes if required. For cover disinfection machine wash at 60ºC using a disinfectant detergent. Do not use softener, chlorine nor whitening concentrates.

WARNING

Max. sliding weight
300 Kg



	Patients <120Kg	Patients 120-200Kg	Patients 200-300Kg
A	20 cm	--	0 cm
B	20 cm	5 cm	0 cm

PACKAGING



- 1 - Medi-roller PRO board
- 1 - Medi-roller PRO cover
- 1 - Medi-roller PRO disposable cover

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE