MEDIROLLER

ENGLISH



Replacement cover for Medi-roller Standard boards *TMR6200 and TMR6220.*

Together with the board, it allows for the lateral movement of patients from one surface to another safely and with minimal effort.

Made of very low friction silicone coated polyamide fabric.

MEDI-ROLLER STANDARD COVER

SPECIFICATIONS

PRODUCT FAMILY: Medi-roller

PRODUCT NAME: Medi-roller Standard cover

INTENDED USE: Replacement sleeve that in combination with the Medi-Roller board, allows patients to be transferred laterally from one surface to another

BASIC UDI-DI: 843603535TMRNE

SIZE CHART

Ref.	Model	Measures (cm)	Fits
TMR6250	Large	180 x 50	TMR6200
TMR6252	Short	90 x 50	TMR6220

COLOUR

Orange

COMPOSITION

Polyamide
LATEX FREE

DISINFECTION RECOMMENDATIONS

- For infection prevention and control use the disposable protective cover (ref. TMR6270).
- For cover disinfection machine wash at 60ºC using a disinfectant detergent. Do not use softener, chlorine nor whitening concentrates.

LAUNDRY INSTRUCTIONS



- Machine wash max. 60ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Medi-roller Standard - cover

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE