MEDIROLLER Pro

ENGLISH



Replacement cover for Medi-roller PRO boards *TMR6205 y TMR6225.*

Together with the board, it allows for the lateral movement of patients from one surface to another safely and with minimal effort.

Made of very low friction polyamide fabric (Ripstop).

MEDI-ROLLER PRO COVER

SPECIFICATIONS

PRODUCT FAMILY: Medi-roller

PRODUCT NAME: Medi-roller PRO cover

INTENDED USE: Replacement sleeve that in combination with the Medi-Roller PRO board, allows patients to be transferred laterally from one surface to another

BASIC UDI-DI: 843603535TMRNE

SIZE CHART

Ref.	Model	Measures (cm)	Fits
TMR6255	Large	165 x 50	TMR6205
TMR6257	Short	80 x 50	TMR6225

COLOUR

- Orange

COMPOSITION

Polyamide
LATEX FREE

DISINFECTION RECOMMENDATIONS

- For infection prevention and control use the disposable protective cover (ref. TMR6275).
- For cover disinfection machine wash at 60ºC using a disinfectant detergent. Do not use softener, chlorine nor whitening concentrates.

LAUNDRY INSTRUCTIONS



- Machine wash max. 60ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Medi-roller PRO - cover

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE