

MEDIROLLER

ENGLISH

TMR6270



TMR6275



***Disposable cover for Medi-roller boards**, Keeps medi-roller boards protected from infections and liquids or body fluids.*

Made of polyethylene plastic supplied in roll of 100 covers precut at 180cm each.

2 different models for fitting Medi-roller Standard or Pro models.

DISPOSABLE COVER

SPECIFICATIONS

PRODUCT FAMILY: Medi-roller

PRODUCT NAME: Disposable cover

INTENDED USE: Protection for Medi-roller boards from potential infections avoiding contact of liquids and body fluids with the board

BASIC UDI-DI: 843603535TMRNE

SIZE CHART

Ref.	Measures (cm)	Fits
TMR6270	180 x 55	TMR6200 - TMR6220
TMR6275	180 x 52	TMR6205 - TMR6225

COLOUR

White

COMPOSITION

Polyethylene
LATEX FREE

LAUNDRY INSTRUCTIONS

- Product for single use
- Do not wash

PACKAGING



- 1 roll - 100 pre-cut disposable covers

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE