# MEDIROLLER

ENGLISH



***Dispenser for disposable covers for Medi-Roller boards,*** *designed for easy handling and cutting of covers. Keeps your Medi-Roller products neatly organized.*

*Stainless steel dispenser with a galvanized steel tubular bar.*

*Comes with a stainless steel wall hook for hanging and storing any Medi-Roller board. The hook is provided separately, allowing you to place it in the most convenient wall location for your needs, above, below, or on the sides of the dispenser.*

## *DISPOSABLE COVER DISPENSER*

## *SPECIFICATIONS*

**PRODUCT FAMILY:** Medi-roller

**PRODUCT NAME:** Disposable cover dispenser

**INTENDED USE:** To hold and dispense covers and medi-roller board. To ease changing cover rolls.

**BASIC UDI-DI:** 843603535TMRNE

### SIZE CHART

| Ref. | Measures (cm) | Fits |
|---|---|---|
| TMR6280 | 600 x 200 x 110 | TMR6270 - TMR6275<br>Mounting hook fits all Medi-roller boards |

### COLOUR

Silver

### COMPOSITION

Dispenser & hook: Stainless Steel
Bar: Galvanized steel

### LAUNDRY INSTRUCTIONS

- Clean with damp cloth and mild soap or surface disinfectant detergent.
- Not recomended to use bleach to desinfect.

### PACKAGING




- 1 Disposable cover dispenser
- 1 Mounting hook
- Wall fixing fittings









MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE