mobi-tools Slide

ENGLISH



TMS6300



TMS6335

***Mobi-tools Padded Tubular Sheet**. Designed to be placed under a bedridden patient to assist the caregiver in repositioning tasks.*

Due to its sliding interior, it also allows the patient to make postural changes by themselves with less effort. The soft padded material provides comfort and can be left in place at all times.

PADDED TUBULAR SLIDING SHEET

FEATURES



Ultra sliding orange fabric



Soft and comfortable padded tubular in wavy stitching finish



5 sizes available

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Slide

PRODUCT NAME: Padded tubular sliding sheet

INTENDED USE: Facilitate repositioning tasks for the caregiver in bedridden patients

BASIC UDI-DI: 843603535TMSNG

SIZE CHART

Ref.	Measures (cm)
TMS6300	25 x 60
TMS6310	50 x 60
TMS6325	100 x 85
TMS6335	140 x 85
TMS6345	190 x 85

COLOUR

- Dark Blue
- Orange

COMPOSITION

Exterior: Polyamide
Padded interior: Polyester
Interior: Polyamide
LATEX FREE

LAUNDRY INSTRUCTIONS




- Machine wash 40ºC
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

WARNING

- Maximum weight allowed: Up to200 Kg.

PACKAGING

- 1 - Mobi-tools Slide - Padded tubular sliding sheet



MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE