mobi-tools Slide

ENGLISH



Mobi-tools tubular one-way *for postural correction of seated patients (wheelchairs, armchairs, car seats, etc.) The unidirectional sliding material prevents the patient from sliding forward on the seat and allows the caregiver to reposition the patient easily and smoothly.*

Made of padded tubular sandwich with a soft touch fabric and anti-slip fabric bands sewn on back for an extra grip. Ultra sliding and unidirectional fabric interior keeps movement in a single direction and blocks the other.

TUBULAR ONE-WAY

FEATURES



Ultra sliding orange fabric



Soft and comfortable padded tubular in wavy stitching finish



Anti-slip bands on back



Interior unidirectional fabric

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Slide

PRODUCT NAME: Tubular one-way

INTENDED USE: For postular correction of seated patients

BASIC UDI-DI: 843603535TMSNG

SIZE CHART

Ref.	Measures (cm)
TMS6360	50 x 45

COLOUR

- Blue/Grey exterior
- Orange/Violet interior

COMPOSITION

Exterior: Polyamide
Padded interior: Polyester
Interior: 80% Polyamide - 20% Cotton
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Tuble dry low
- Do not iron
- Do not dry clean

WARNING

- Maximum weight allowed: Up to 200 KG.

PACKAGING



- 1 - Mobi-tools Slide - Tubular one-way

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE