mobi-tools Grip

ENGLISH



Mobi-tools mobilisation band *provides the therapist or caregiver with extra extension arms. Allows controlling the patient trunk in mobilisation and transfer tasks, providing a sense of security to the patient and avoiding overexertion injuries in the caregiver.*

Made with an 8mm polyester foam textile sandwich, with one side of polyester canvas and an interior of easy-to-clean PVC fabric. High-strength polyester straps.

MOBILISATION BAND

FEATURES



2 grip positions on each side for transfer and mobilization tasks



Padded belt with an inner PVC fabric that is easy to clean

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Grip

PRODUCT NAME: Mobilisation band

INTENDED USE: Ergonomic and safe patient control during mobilization and transfer tasks

Basic UDI-DI: 843603535TMTNJ

SIZE CHART

Ref.	Size	Measures (cm)	
		Length	Width
TMT6500	One size	87	16

COLOUR

- Orange
- Dark blue

COMPOSITION

Exterior: 70%PES - 30%PVC
Interior: Polyester
Straps: Polyester
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Mobi-tools Grip - Mobilisation belt

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE