mobi-tools Grip

ENGLISH



Mobi-tools semirigid mobilisation band *provides the therapist or caregiver with extra extension arms. Allows controlling the patient trunk in mobilisation and transfer tasks, providing a sense of security to the patient and avoiding overexertion injuries in the caregiver.*

Incorporates a plate inside to facilitate its placement between the patient and the bed or chair, preventing it from bending. Made with resistant materials and double handles at the ends with reinforced seams ensure safe, effective and durable use.

SEMIRIGID MOBILISATION BAND

FEATURES



2 Padded handles straps on each side



Comfortable padded soft to the touch interior fabric



Rigid body to facilitate the placement between patient and bed

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Grip

PRODUCT NAME: Semirigid mobilisation band

INTENDED USE: Ergonomic and safe patient control during mobilization and transfer tasks

Basic UDI-DI: 843603535TMTNJ

SIZE CHART

Ref.	Size	Measures (cm)	
		Length	Width
TMT6505	M	50	20
TMT6506	L	60	20

COLOUR

- Orange
- Dark blue

COMPOSITION

Exterior: 70%PES - 30%PVC
Interior: Polyester
Straps: Polyester
Reinforcement:
3mm PEHD plate / 2mm EVA Foam
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Mobi-tools Grip - Semirigid mobilisation belt

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE