mobi-tools Grip

ENGLISH



Mobi-tools mobilization belt. *It can be placed on the patient or the caregiver, providing an extra support element for the mobilization and transfer of patients. In this way, a greater distribution of effort is achieved, avoiding lumbar overload of the caregiver and increasing the patient's feeling of security. It adapts to the waist circumference via a sliding buckle. Ergonomically placed handles are provided for transfer and mobilization tasks, as well as dorsal handles for walking assistance.*

Made with an 8mm polyester foam textile sandwich, featuring one face of polyester canvas and an interior of PVC fabric for easy cleaning. 50mm high-resistance polyester straps. Acetate (POM) buckles.

MOBILISATION BELT

FEATURES



Grip handles on each side, ergonomically arranged for transfer and mobilization tasks



Padded belt with PVC inner fabric for easy sanitation



Horizontal handles on back for deambulation

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Grip

PRODUCT NAME: Mobilisation belt

INTENDED USE: Provides extra support for mobilization and patient transfer tasks

Basic UDI-DI: 843603535TMTNJ

SIZE CHART

Ref.	Size	Measures (cm)	
		Contour (Min/Max)	Width
TMT6511	S	60 / 90	10
TMT6512	M	80 / 120	14
TMT6513	L	90 /140	14

COLOUR

- Orange
- Dark blue

COMPOSITION

Exterior: 70%PES - 30%PVC
Interior: Polyester
Straps: Polyester
Acetate (POM) buckles
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Mobi-tools Grip - Mobilisation belt

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE