mobi-tools Grip

ENGLISH



Mobi-tools lower legs band *it provides an extra grip, support and anchoring point at the bottom of the person's legs. Very useful in transfer and mobilisation tasks that require the patient's legs to be moved, for instance, when helping them get out or into the bed, or when blocking their legs to stand up or sit down.*

Made with an 8mm polyester foam textile sandwich, featuring one face of polyester canvas and an interior of PVC fabric for easy cleaning. High-strength polyester straps. Easy and quick adjustment using a hook closure over resin-coated velour fabric.



LOWER LEGS BAND

FEATURES



Side handles for positioning



Padded with an inner PVC fabric that is easy to clean



Quick and easy hook/loop adjustment

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Grip

PRODUCT NAME: Lower legs band

INTENDED USE: Provides extra support for mobilization and patient transfer tasks

BASIC UDI-DI: 843603535TMTNJ

SIZE CHART

Ref.	Measures (cm)
TMT6520	110 x 24

COLOUR

- Orange
- Dark blue

COMPOSITION

Exterior: 70%PES - 30%PVC
Interior: Polyester
Straps: Polyester
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 - Mobi-tools Grip - Lower legs band

MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE