mobi-tools Grip

ENGLISH



Mobi-tools Assistive Ladder. *It consists of an aid device for the incorporation of bedridden or wheelchair-bound patients.*

It has 6 separate bands in different colors to make identification easy. Each handle incorporates a rubber grip for comfortable pulling.

Made with highly durable and strong polyester straps, recommended for patients up to 200 kg.

PULL & GRIP BAND

FEATURES



High strength and durability polyester straps



Easy to install to the bed frame



Soft rubber grip



Patient weight up to 200 Kg.

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Grip

PRODUCT NAME: Pull & Grip band

INTENDED USE: Facilitate the incorporation of the bedridden or wheelchair patient

Basic UDI-DI: 843603535TMTNJ

SIZE CHART

Ref.	Measures (cm)
TMT6540	200 x 20

COLOUR

- Blue
- Green
- Red
- Yellow

COMPOSITION

Polyester
Rubber grips
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 70ºC
- Do not bleach
- Tumble dry low
- Iron low
- Do not dry clean

PACKAGING



- 1 Mobi-tools Grip - Pull & Grip band

WARNING

- Maximum weight allowed: Up to200 KG.

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE