
mobi-tools Turn

ENGLISH



***Mobi-tools rigid Turntable**. Allows smooth, safe, and controlled rotation of the patient during transfer movements requiring rotation.*

It is covered with anti-slip material on both sides to adhere to all types of surfaces and provide more stability and safety for the patient. It incorporates a puller handle to facilitate gripping the disc from the floor. Height of only 15 mm.

RIGID TURNTABLE

FEATURES



Smooth & effortless rotation



Top & Bottom anti-slip surface



Patient weight up to 135 Kg.

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Turn

PRODUCT NAME: Rigid turntable

INTENDED USE: Facilitate the rotation of bedridden or wheelchair patient

Basic UDI-DI: 843603535TMTNJ

SIZE CHART

Ref.	Measures (mm)		Weight (Kg)
	Diameter	Height	
TMT6550	Ø 415	14,5	1,68

COLOUR

- White
- Grey anti-slip

COMPOSITION

Lower disc: ABS
Upper disc: PVC
Anti-slip: 50% PU - 50% PES
LATEX FREE

LAUNDRY INSTRUCTIONS

- Hand wash using a damp cloth or brush in an aqueous solution and neutral detergent.
- Use a surface disinfectant or sanitizing disinfectant wipes, if required.

WARNING

- Maximum weight: 135 Kg.

PACKAGING



- 1 - Mobi-tools Turn - Rigid turntable

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE