mobi-tools Turn

ENGLISH



Mobi-tools swivel cushion *facilitates the rotation of patients or dependent person on the bed, chair or car seat.*

Manufactured in a combination of 8mm foam bonded with a soft touch polyester fabric. Inside it incorporates an ultra-sliding and very low friction polyamide technical material which facilitates rotation. Bottom face made in anti-slip coated fabric. Two sewn polyamide pieces work as an axis and as joint.

SWIVEL CUSHION

FEATURES



Smooth & effortless rotation



Bottom anti-slip surface



Ultra sliding fabric inside

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Turn

PRODUCT NAME: Swivel cushion

INTENDED USE: Facilitate the rotation of bedridden, wheelchair patients or car seat patient

Basic UDI-DI: 843603535TMTNJ

SIZE CHART

Ref.	Measures (mm)	
	Diameter	Height
TMT6555	Ø 400	25

COLOUR

- Blue/Grey (exterior)
- Orange (interior)

COMPOSITION

Interior: Polyamide
Exterior: Polyester
Turning plate: Polyamide
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash max. 40ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

WARNING

- Maximum weight: 135 Kg.

PACKAGING



- 1 - Mobi-tools Turn - Swivel cushion

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE