mobi-tools Translide

ENGLISH



Mobi-tools tubular sliding sheets *are a useful aid to carry out tasks of repositioning patients in bed and as complements for performing transfers. Made of high tenacity silicone coated nylon fabric, that gives it super-slippery attributes and helps caregivers to perform their jobs with greater comfort and safety, preventing strain injuries due to low friction. Sewn loop strap for wall hanging.*







TUBULAR SLIDING SHEET

SPECIFICATIONS


PRODUCT FAMILY: Mobi-tools Translide

PRODUCT NAME: Tubular sliding sheet

INTENDED USE: Aid to carry out tasks of repositioning patients in bed and as complements for performing transfers

BASIC UDI-DI: 843603535TTEN9

SIZE CHART



Ref.	Measures (cm)	
	A	B
TTE6100	60	41
TTE6102	72	71
TTE6104	122	71
TTE6106	122	100
TTE6108	200	71

COLOUR

Red

COMPOSITION

Polyamide
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

WARNING

- Maximum weight: 150 Kg.


HYGIENIZATION / DISINFECTION

Machine wash at a temperature of 60ºC using a low-temperature disinfection detergent

Superficially with surface disinfectants or disinfectant wipes

PACKAGING

- 1 - Mobi-tools Translide - Tubular sliding sheet



MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE