mobi-tools Translide

ENGLISH



TTE6110-6112



TTE6122

Mobi-tools sliding sheet *is a convenient solution for transfer and mobilization of bedridden patients. For safety transfers it is recommended their combined use with Medi-glide boards TMG6420 or TMG6430.*

Made of ultra-sliding high-tenacity low friction technical fabric that helps caregivers to perform their jobs preventing strain injuries.

TTE6122 *- Includes handles on the full length and the top side of the sheet (3 sides) providing the healthcare staff more control of the patient and maintain the correct posture during the movement.*

SLIDING SHEET

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Translide

PRODUCT NAME: Sliding sheet

INTENDED USE: Aid to carry out tasks of repositioning patients in bed and as complements for performing transfers

BASIC UDI-DI: 843603535TTEN9

SIZE CHART

Ref.	Measures (cm)		Handles
	Length	Width	
TTE6110	200	70	-
TTE6112	200	140	-
TTE6122	200	90	3 sides

COLOUR

- Red

COMPOSITION

Sheet: Polyamide
Handles (TTE6122): Polypropylene
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

WARNING

- Maximum weight: 150 Kg.

HYGIENIZATION / DISINFECTION

Machine wash at a temperature of 60ºC using a low-temperature disinfection detergent

Superficially with surface disinfectants or disinfectant wipes

PACKAGING

- 1 - Mobi-tools Translide - Sliding sheet



MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE