mobi-tools Translide

ENGLISH



Mobi-tools bariatric transfer sheet *is a convenient solution for transfer and mobilization of bedridden patients. For safety transfers it is recommended their combined use with Medi-glide boards TMG6420 or TMG6430.*

Made of ultra-sliding high-tenacity low friction technical fabric with 8 handles, 4 on each side that helps caregivers to perform their jobs preventing strain injuries.

BARIATRIC TRANSFER SHEET

FEATURES

X-RAY / MRI — **X-ray translucent and MRI safe**

Silicone coated polyamide fabric

4 handles on each side

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Translide

PRODUCT NAME: Bariatric transfer sheet

INTENDED USE: Aid to carry out tasks of repositioning patients in bed and as complements for performing transfers

BASIC UDI-DI: 843603535TTEN9

SIZE CHART

Ref.	Measures (cm)		Handles
	Length	Width	
TTE6142	200	90	8

COLOUR

- Blue
- White handles

COMPOSITION

Sheet: Polyamide
Handles: Polypropylene
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 70ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

WARNING

- Maximum weight: 180 Kg.

HYGIENIZATION / DISINFECTION

Machine wash at a temperature of 70ºC using a low-temperature disinfection detergent

Superficially with surface disinfectants or disinfectant wipes

PACKAGING

- 1- Mobi-tools Translide - Bariatric transfer sheet



MEDICARE SYSTEM
MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE