mobi-tools Translide

ENGLISH



Mobi-tools sliding mittens *are an effective aid to reduce friction when handling the patient during repositioning or mobilisation tasks. It also allows to accommodate the patient in a sitting position and eliminate wrinkles in the sheet that may cause discomfort. Made of ultra sliding silicone/PU coated Nylon fabric .*





SLIDING MITTENS

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Translide

PRODUCT NAME: Sliding mittens

INTENDED USE: Aid to carry out tasks of repositioning patients in bed and as complements for performing transfers

BASIC UDI-DI: 843603535TTEN9

SIZE CHART

Ref.	Measures (cm)	
	Width	Length
TTE6150	20	50

COLOUR

Blue

COMPOSITION


Polyamide
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 70ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING



- 1 Pair - Mobi-tools Translide - Sliding mittens

HYGIENIZATION / DISINFECTION

Machine wash at a temperature of 70ºC using a low-temperature disinfection detergent

Superficially with surface disinfectants or disinfectant wipes

MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE