mobi-tools Translide

ENGLISH



Mobi-tools disposable transfer sheet *is a convenient solution for transfer and mobilization of bedridden patients.*

For safety transfers it is recommended their combined use with Medi-glide boards TMG6420 or TMG6430.

Made of water repellent polypropylene non-woven fabric with 5 handles on each side that helps caregivers to perform their jobs preventing strain injuries.

DISPOSABLE TRANSFER SHEET

SPECIFICATIONS

PRODUCT FAMILY: Mobi-tools Translide

PRODUCT NAME: Disposable transfer sheet

INTENDED USE: for transfer and mobilization of bedridden patients

BASIC UDI-DI: 843603535TTEN9

SIZE CHART

Ref.	Measures (cm)		Handles	Weight (gr/m²)
	Length	Width		
TTE6160	194	85	10	95

COLOUR

- Dark blue

COMPOSITION

Polypropylene
LATEX FREE

WARNING

- Max. Weight: 150 Kg.


LAUNDRY INSTRUCTIONS



- Do not wash
- Disposable Product

PACKAGING

- 4 - Mobi-tools Translide - Disposable transfer sheet



MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE