
mobi-tools Translide

ENGLISH



Mobi-tools tubular sliding sheets for wet environment *are a useful aid to carry out tasks of repositioning patients in bed and as complements for performing transfers.*

For safety transfers it is recommended their combined use with Medi-glide boards TMG6420 or TMG6430.

Made of high tenacity ripstop silicone coated nylon fabric on both sides, that gives it super-slippery attributes and helps caregivers to perform their jobs with greater comfort and safety, preventing strain injuries due to low friction. Fabric specially designed for wet environments, non-porous and sliding in wet conditions.

TUBULAR SLIDING SHEET (WET ENVIRONMENTS)

SPECIFICATIONS


PRODUCT FAMILY: Mobi-tools Translide

PRODUCT NAME: Tubular sliding sheet (wet environments)

INTENDED USE: aid to carry out tasks of repositioning patients in bed and as complements for performing transfers in wet environments

BASIC UDI-DI: 843603535TTEN9

SIZE CHART



Ref.	Measures (cm)	
	A	B
TTH6020	100	66
TTH6040	150	66

COLOUR

Blue

COMPOSITION

Polyamide
LATEX FREE

WARNING

- Maximum weight: 150 Kg.

LAUNDRY INSTRUCTIONS



- Machine wash 60ºC
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING

- 1 - Mobi-tools Translide - Tubular sliding sheet (Wet)



MEDICARE SYSTEM

MEDI CARE SYSTEM S.L.U.
C/ Abat Escarré 15 - 08302 Mataró SPAIN
Tel.: +34 937 570 048 - Fax: +34 937 414 068
E-mail: info@medisystemsl.com · www.medicaresystem.es

CE Declaration according to Regulation (EU)745/2017 of the European Parliament & Council

*This data sheet could be changed without prior notice

CE